




**Status of BL13-XALOC
Macromolecular Crystallography
beamline**



Jordi Juanhuix

Beamline concept : In-vacuum Undulator IVU21.6 + Channel-cut Si(111) monochromator + KB focusing mirrors








Experimental
Hutch

MD2M
CATS
PILATUS 6M

Beam size be expanded by unbending the mirrors and defocusing the beam

- From 50×5.7 to 250×250 μm
- Increasingly popular for beamtimes with a variety of projects




Beam Shaping - smaller sizes

The MD2 penta-aperture consisting in 5 round pinholes is used to reduce the beam size


- Variety of smaller beam sizes for small crystals
- More accurate raster scans possible

BUT:



- Flux is reduced linearly with size
- Data have slightly but consistently better quality (ISAs) when defocusing the beam
- Seeking for beam vibrations...



V profile at sample




H profile at sample





Beam vibrations

- Beam center found by Gaussian fitting of images at focal point recorded at 480 Hz
- High speed camera coupled with a 20 μm -thick YAG:Ce screen
- Vibrations with the beamline at rest (no excitations)


Horizontal oscillation, $\sigma = 0.24 \mu\text{m}$




Vertical oscillation, $\sigma = 0.35 \mu\text{m}$




- The origin of the resonances can be found by hitting the vacuum vessels of the critical elements
- All resonances have frequencies higher than the maximum detector frame rate (12 Hz)




Mono eigenfrequencies



VFM eigenfrequencies



HFM eigenfrequencies






- Horizontal beam vibrations are dominated by a 48 Hz coming from HFM
- Vertical beam vibrations are dominated by a series frequencies at ~75 Hz and one at 100 H from Mono and VFM
- Still no conclusions on the effect on the quality of data.





Spectrum of the beam motion at detector position (no-excitation): Amplitudes

ID	frequency interval (Hz)	frequency (Hz)	x-amplitude (um rms)	ratio to total amplitude (%)	y-amplitude (um rms)	ratio to total amplitude (%)		Resonance
1	0.00-15.00	6.53	0.034	2.08	0.06	3.05		
2	15.00-22.00	19.16	0.01	0.17	0.029	0.71		
3	22.00-35.00	26.13	0.011	0.22	0.026	0.55		
4	35.00-41.00	39.42	0.011	0.21	0.014	0.15		
5	41.00-45.00	43.2	0.014	0.35	0.019	0.31		VFM
6	45.00-49.60	48.71	0.164	47.66	0.017	0.23	H	HFM
7	49.60-52.00	50.04	0.038	2.55	0.021	0.38		
8	52.00-67.00	58.53	0.042	3.14	0.113	10.73		
9	67.00-72.00	69.6	0.048	4.05	0.092	7.15		
10	72.00-78.00	73.47	0.121	25.83	0.131	14.37	Both	HFM, VFM, MONO
11	78.00-88.00	78.8	0.036	2.25	0.102	8.78		
12	88.00-98.00	95.82	0.016	0.43	0.058	2.82		Mono
13	98.00-106.00	100.36	0.036	2.25	0.165	22.9	V	
14	106.00-110.00	108.93	0.009	0.15	0.043	1.57		VFM
15	110.00-116.00	112.93	0.01	0.18	0.054	2.4		
16	116.00-124.00	116.98	0.013	0.31	0.09	6.82		
17	124.00-129.00	126.27	0.009	0.13	0.038	1.18		
18	129.00-142.00	137.42	0.017	0.48	0.035	1.04		
19	142.00-145.00	143.2	0.008	0.11	0.02	0.32		
20	145.00-149.00	148.27	0.009	0.15	0.053	2.33		
21	149.00-150.50	150.04	0.037	2.4	0.029	0.71		
22	150.50-154.00	151.24	0.011	0.21	0.049	2.05		
23	154.00-160.00	156.13	0.012	0.26	0.038	1.24		Mono
24	160.00-168.00	164.18	0.021	0.75	0.056	2.66		
25	168.00-200.00	169.6	0.046	3.74	0.082	5.7		

- Plates are mounted on the robot, using the plates tool
- The new XYZ moveable beam stop required for plates is already installed
- The cryostream is swung out remotely
- A special widget is used for collection. Integration in the standard collection widget is foreseen before Summer 2016
- A TRIS crystal was positioned in the beam using manual operation of the robot
- An still image was recorded using an exposure of 0.5 s
- Now able to scan plates, but not ready for users




- **EDNA** implemented and integrated in the experiment pipeline
- **Automatic data processing**
 - Using the cluster
 - AutoPROC (Global Phasing Ltd) implemented
 - Implementing Phaser in the pipeline
- Diffraction **data queues** ready.
 - **Inverse beam** (wedged) anomalous data collection
 - Many others to be implemented
- Multiple crystal position data collection: **crystal cartography**
- **Reorientation of the crystal** after indexing with EDNA
- **Beamline diagnostic** tool
- **sftp service for users**



- NX based (Server ready)
- User needs to install NX Enterprise client for Windows/Linux/Mac (exe file)
- Connect to:
Host: `remotensex.cells.es`
Port: 443
- Currently in beta testing
- Still to be done: script giving user authentication

WORKFLOW REMOTE ACCESS XALOC





Mechanical

- **Installation of a (de)humidifier device**
 - Control system required
- **Mechanized crystal annealing**
- **New viewing camera**



Controls

- **Installation of MXCuBE**
- Implementation of (pending MxCuBE??):
 - **Helical** data collection
 - **MAD** data collection **pipeline**
 - **Multiple kappa** angle data collection
 - Automated crystal centering
- **Automatic beam alignment** to be started when the fluorescence screen can be automatically placed at sample position
- **Automatic loop centering** : waiting for MXCuBE?



Engineering

Nahikari González
Carles Colldelram
Robert Oliete
Ricardo Valcárcel

Mech. Engineer
Mech. Engineer
Mech. Technician
Mech. Technician

Computing

Jordi Andreu
José Ávila
Xavi Fariña
Jorge Villanueva

Software Engineer
Electronics Engineer
Electronics Tech.
EPS engineer

Experiments

Fernando Gil
Roeland Boer
Jordi Juanhuix
Josep Nicolás
Igor Šics

BL scientist
BL scientist
BL responsible
Mirrors, Metrology
Diagnostics

Inma Hernández

Administration

